

Supplementary Appendix

This appendix has been provided by the authors to give readers additional information about their work.

Supplement to: McMahon MA, Jilek BL, Brennan TP, et al. The HBV drug entecavir — effects on HIV-1 replication and resistance. *N Engl J Med* 2007;356:2614-21.

Supplemental Methods

Real Time PCR Primers

The following set of primers and probes were used for Real Time PCR amplification:

5' early, GCTAACTAGGGAACCCACTGCTT; 3' early,
CAACAGACGGGCACACACTACT; early probe,
AGCCTCAATAAAGCTTGCCTTGAGTGCTTC;

5' intermediate, GGTGCGAGAGCGTCGGTATTAAG; 3' intermediate,
AGCTCCCTGCTTGCCCATA; intermediate probe,
AAAATTCGGTTAAGGCCAGGGGAAAGAA;

5' late, GAGTCCTGCGTCGAGAGATC, 3' late, TGTGTGCCCGTCTGTTGTGT; late
probe, CAGTGGCGCCCGAACAGGGA.